

**Product Name :**  
Cryo Vial Internal Thread

**Product Code :**  
CHEMISTRYSUK-AS03LAB009R



**Description :**

Cryo Vial Internal Thread

**Technical Specification :**

The cap fitted with an O-ring screw into the vials and makes it the leakproof Cryo Vial.

Adding to the existing range of Cryo Vials is this which has internal thread in it.

These vials are provided with white, prominent graduations and a long marking area for sample identification.

**Capacity:** 2.0ml



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265