

Product Name :
Cryo Vial

Product Code :
CHEMISTRYSUK-AS03LAB010R



Description :

Cryo Vial

Technical Specification :

The cap of these vials is provided with a seal ring inside to make it completely leakproof.

Cryo vials are normally used for storage of Biological material, human & animal cells at temperature as low as -190°C.

These sterilized tubes are provided with white, prominent graduations and a long marking area for sample identification.

The longitudinal grooves in the base give a non twisting feature to these vials when placed in a cryo rack.

Capacity: 4.5ml, 1.8ml, 1.0ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265