

**Product Name :**  
Friability Test Apparatus

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB001A



**Description :**

Friability Test Apparatus

**Technical Specification :**

The Friability Test Apparatus is designed as per IP/USP standards. The unit is equipped with two transparent acrylic drums which rotate at a speed of 25 . 1 rpm. One transparent acrylic drums rotated by superior quality constant speed sturdy motor, mounted in painted sheet metal cabinet with chemical resistant powder coating. The apparatus is designed to determine the abrasion rate & impact hardness of tablet. Friability Test Apparatus is a precision instrument for measuring obression strength of the tablets.

**Specification:-**

Speed: 25 + 1 RPM

Motor: 1/8 hp

Power: AC/ 230 Volts 50 Hz Single Phase



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265