

**Product Name :**  
Dissolution Six Test Apparatus

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB002A



**Description :**

Dissolution Six Test Apparatus

**Technical Specification :**

The round bottom glass jars are placed equidistant inside the motors. The motors are electronically controlled to speed precisely within +/- 1 RPM. Dissolution Six Test Apparatus is designed to carry out test of tablet and capsule as per USP Standard. These tests are carried out by quality control and R & D departments in pharmaceutical industry. It consists of constant temperature water bath made of transparent acrylic. Temp. is controlled by a fully solid-state controller within +/- 0.5 deg C. The RPM is controlled digitally and is indicated on a digital display. Motor assembly can be lifted with an auto lifting mechanism for easy operation. The basket assemblies and stirring paddles are made of 316-grade stainless steel. It is a microprocessor based digital pre-set timer with audio alarm.

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

