

Product Name :
Conical Percolator

Product Code :
PHARMACEUTICAL SUK-
AS03LAB004A



Description :

Conical Percolator

Technical Specification :

These percolators are provided with lid drain cock which is supported on a stainless steel sieve on 3-legged stand. Conical Percolator is made up of stainless steel. Our percolator is widely used in laboratory experiments which are conducted for research as well as in teaching purposes in institutional and research centers.

Features :-

Accurate Dimension.

Anti Corrosive Body.

Long Service Life.

Simple and Safe in Use.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265