

Product Name :
Tincture Press

Product Code :
PHARMACEUTICAL SUK-
AS03LAB005A



Description :

Tincture Press

Technical Specification :

Press is built with stainless steel 304 material. Easy to operate and low maintenance. Capacity of the press ranges from 50-200 liter. Ensure longer working life. Provides high quality tincture.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265