

**Product Name :**  
Hand Operated Tablet Machine

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB009A



**Description :**

Hand Operated Tablet Machine

**Technical Specification :**

The machine is supplied without dies / punches. Hand Operated Tablet Making Machine is driven hand operated by circulating wheel. The compression pressure on the tablets, is controlled by adjusting a lever, the powder filled in cavity is controlled automatically by a moving a feeder provided. Complete with Hopper, feeder, Tablets slider etc.

**Features :-**

Compressing lever for adjustable harness of tablets.

Adjustable lower punch for adjusting the thickness.

Easily rotating wheel for making tablets.

Stainless steel day and punch size 6, 9, 12 mm.

Electrical model can be provided in same with ½h.p motor.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265