

**Product Name :**  
Tablet Coating Pan

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB010A



**Description :**

Tablet Coating Pan

**Technical Specification :**

With the aid of proficient production engineers and allied professionals, we manufacture these tablet coating pans in compliance with international standards. Tablet Coating Pan Stainless steel is used in manufacturing the pan whereas rest of the machine parts are made out of mild steel.

**Tablet coating pans consist of the following :-**

Thermostatically controlled heaters.

Complete with suitable motor.

0.5 HP hot air blower with suitable starter and indicator.

Variable pulley for controlling of R.P.M of moving pan.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265