

Product Name :
Polishing Pans

Product Code :
PHARMACEUTICAL SUK-
AS03LAB011A



Description :

Polishing Pans

Technical Specification :

The pans can be fitted with internal baffles if required. This helps tumble the product for better coating. Due to their simple design the machines are very easy to operate and maintain. The angle of the pan can be adjusted to give better results and capacities on different products. The pans can also be supplied with hot or cold air blower systems or automatic spray systems.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265