

**Product Name :**  
Tablet Disintegration Test Apparatus

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB014A



**Description :**

Tablet Disintegration Test Apparatus

**Technical Specification :**

The Instrument is also equipped with double channels Temperature system to ensure precise temperature control of the Tank and Beakers. Illuminator for clear observations Of Disintegration process has been provided. An instrument is equipped with the single stepper motor. Five programmable Tests are provided to Store DT time and Temperature. Digital indicator for temperature with optically isolated electronic heater driver with the water bath. Soft membrane keypad and splash proof panel. A four line informative LCD display is provided to observe temperature and set time at a glance. The facility provided to prevent over the shooting of temperature. Facility for self-validation of RPM is incorporated. Automatic Parking of both Baskets after the Test is over. The equipment meets all international standards with respect to Calibration, Validation and IQ, OQ requirement. DT-TIME for Tablets from individual Tube Of Right and Left Basket can be monitored. Centronics Printer interface has been provided for recording parameters like Temperature, Batch No, Product ID etc.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265