

**Product Name :**  
Collapsible Tube Filling and Sealing Machine

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB018A



**Description :**

Collapsible Tube Filling and Sealing Machine

**Technical Specification :**

Built on cast iron body, hand lever operated. Stroke type, barrel and piston arrangement with stroke setting for various capacity of filling of tubes (from 5 gms to 50 gms.) detachable nozzles provided for each separate capacity. All contact parts are of stainless steel also stainless steel holder.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265