

Product Name :
Collapsible Tube Filling and Sealing Machine

Product Code :
PHARMACEUTICAL SUK-
AS03LAB018A



Description :

Collapsible Tube Filling and Sealing Machine

Technical Specification :

Built on cast iron body, hand lever operated. Stroke type, barrel and piston arrangement with stroke setting for various capacity of filling of tubes (from 5 gms to 50 gms.) detachable nozzles provided for each separate capacity. All contact parts are of stainless steel also stainless steel holder.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265