

**Product Name :**  
Capsule Filling Machine

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB019A



**Description :**

Capsule Filling Machine

**Technical Specification :**

Handle locking arrangement makes the operator to lift the loading plate without fear of damaging the capsules. Main body made of heavy cast iron and all contact parts made of brass richly plated or stainless steel. A compact hand operated machine suitable to fill two-part hard gelatine capsule .

Plate of stainless steel for any one no. Is supplied with machine.

Pressing is achieved by cam system which also saves time to press the plunger as compared to screw system.

For filling 100 capsules.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265