

Product Name :
Distillation Unit For Distilled Water

Product Code :
PHARMACEUTICAL SUK-
AS03LAB021A



Description :

Distillation Unit For Distilled Water

Technical Specification :

It is easy to use, safe to operate. Distilled Water designed for reliable continuous operation incorporating a host of features unmatched by comparable stills. Abundant Output. The distilled output is cool ensured by a high-efficiency condenser. The still provides, 4 Liter/hr of distilled water produced through a power input of 3KW by a chromium-plated heater housed in a horizontal Glass Boiler. The temperature of distillate 25o C – 40o C tested for an ambient room temperature of 30o C.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com