

**Product Name :**  
Glass Distillation

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB023A



**Description :**

Glass Distillation

**Technical Specification :**

These apparatus are used for single, double and triple distillation in various capacities. The Automatic Glass Distillation Apparatus made accessible by us is available with metal stand, electric plug and board. An automatic cut off device is also attached with glass distillation apparatus. In addition, we are one of identified as one of the eminent Glass Apparatus Exporters from the country.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265