

Product Name :
Ampoule Washing Device

Product Code :
PHARMACEUTICAL SUK-
AS03LAB024A



Description :

Ampoule Washing Device

Technical Specification :

Spring loaded ampoule holding arrangement for various sizes of ampoules and a stainless steel bottom tray on the stand. Jet needle type, consists of 6 stainless steel needle, jet fixed on a stainless steel pipe of 12.5 mm. Dia. With side tap and serrations for fixing rubber tube with tap water.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265