

Product Name :
Ampoule Filling Sealing Device

Product Code :
PHARMACEUTICAL SUK-
AS03LAB025A



Description :

Ampoule Filling Sealing Device

Technical Specification :

Electrically operated with piston of syringe having adjustment for quantity of liquid to be filled. Completed with burner for sealing. Possibility of gas flushing before, during and after filling. Ampoule gas centering during gas flushing and filling. Fast change over due to extremely limited parts to be replaced. Possibility to memorize various parameters of different ampoule sizes. Extremely robust, high engineered.

Machine equipped with pre-warming station.

Minimum change over time, little maintenance, easy cleaning. AISI 316L dosing syringes.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265