

**Product Name :**  
Autoclave

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB028A



**Description :**

Autoclave

**Technical Specification :**

**Specification :-**

- Chamber Type: Round Vertical
- Temperature Controller: Digital PID Controller
- Chamber Size: 350 X 630mm

- Capacity: 40 Liter
- Temperature Range: Ambient to 132C
- Operating Pressure: 1.2 kg/cm2 at 121C
- Timer: 99 min. Digital Timer
- Pressure Gauge: Mechanical Type 0 to 3 kg/cm2
- Air Exhaust: Adjustable Handle Valve
- Safety Device: Over Heat Protector, Over Pressure Protector, Alarm
- Power: 3 kw
- Voltage: 220V AC, 60Hz

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com