

**Product Name :**  
Hot Air Sterilizer Chamber

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB029A



**Description :**

Hot Air Sterilizer Chamber

**Technical Specification :**

Exterior fabricated out of thick mild steel sheet duly finished in stove enamel off-white paint or paste mat finish colour combinations. The equipment is suitable to operate on 220 V AC 50 Hz single phase. Inner chamber made of richly anodized aluminum or highly polished stainless steel sheet. The gap between the walls is filled with-special grade glass-wool for heat insulation. Inner chamber is ribbed for placing shelves at convenient levels.

Nichrome wire heating elements -provided on three sides attain quick and uniform heating in range of 50oC to 250oC + 1oC controlled by capillary type thermostat. Provided with two or three removable shelves. Special elegant door closing device has been provided. L-shaped thermometer is built-in type. Control panel is provided with selector switch of high or low rates of power, thermostate control knob and indicators for mains & thermostat.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265