

**Product Name :**  
Laboratory Blender

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB032A



**Description :**

Laboratory Blender

**Technical Specification :**

The container is stainless steel with stainless steel blades and a stainless steel clamp-on lid with rubber gasket and stopper. Laboratory blender features a sealed, precision ball bearing system for a long service life. 3.75 HP motor offers 3 speeds plus pulse. High speed is 20,800 rpm, medium is 18,000 rpm and low is 15,800 rpm. The deeper, wider, full-traction rubber, non-skid feet keep the unit stable. Electrical rating is 120VAC, 50/60 Hz. A 240-volt model is available.

**Features:-**

Easy-to-clean electronic membrane panel

Two handles for safer, easier lifting, pouring and carrying

Removable dishwasher-safe jar pad



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265