

Product Name :
Lab Centrifuge 4 Tubes

Product Code :
PHARMACEUTICAL SUK-
AS03LAB033A



Description :

Lab Centrifuge 4 Tubes

Technical Specification :

This mini centrifuge is great for all types of laboratories, especially clinical and research labs. Widely used in physical, biochemistry, cellular and molecular biology, clinical labs, blood donation centers and much more. The extremely stable with little to no noise. It is built with the highest quality brushless DC motor for the best performance and longest life possible. High-performance maintenance-free brushless motor. LCD display and control of time and speed.

Speed: 300 - 5000 rpm / Timer 15s - 99min.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265