

Product Name :
Refractometer

Product Code :
PHARMACEUTICAL SUK-
AS03LAB037A



Description :

Refractometer

Technical Specification :

The Digital refractometer with a digital readout and LCD display will eliminate uncertainties and is more convenient to use. The Brix concentration, Refractive index or other parameters are always needed in the lab or our life. An optical refractometer is very popular used to measure these parameters. .

Features :-

Clear vision, easy to operate

Factory Calibrated

Economic, best for field use

Available in different ranges starting from zero to 92%

Maintenance free



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265