

Product Name :
PH Meter Digital Table Model

Product Code :
PHARMACEUTICAL SUK-
AS03LAB040A



Description :

PH Meter Digital Table Model

Technical Specification :

These instruments are usually used to measure the pH of liquids, though special probes are sometimes used to measure the pH of semi-solid substances. It is a device used for correctly measuring the pH, which is either the concentration or the activity of hydrogen ions, of an aqueous solution. It usually has a glass electrode plus a calomel reference electrode, or a combination electrode.

Features :-

Easy to use

Special probes

Used in pH measurement



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265