

**Product Name :**  
Haemoglobinometer

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB043A



**Description :**

Haemoglobinometer

**Technical Specification :**

The two side windows are fitted with colour standard strips and middle window is meant for the square tube for colour comparison. The back of the holder is fitted with milky colour screen made of acrylic. In front of the holder there are three windows. The size of all the three windows are Height 3 cms. and Width 0.5 cms.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265