

Product Name :
Student Organ Bath

Product Code :
PHARMACEUTICAL SUK-
AS03LAB045A



Description :

Student Organ Bath

Technical Specification :

The bath size 8"x4"x5" with leak proof dovetailed joints, and rests on Stainless Steel supports. The Student organ bath with thermostatic temperature controlled has been designed mainly for student experiments with smooth muscle etc. Capacity of jar approximately 3/4 liters. Supplied with Uprights, organ & Oxygen tubes with platinum tip, frontal lever from SS capillary tubing.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265