

Product Name :
Sphygmomanometer

Product Code :
PHARMACEUTICAL SUK-
AS03LAB057A



Description :

Sphygmomanometer

Technical Specification :

Adjustable display column for easier analysis of result. Anti-glare display minimizes reflected light. Auscultatory blood pressure measurement. Digital readout for pulse rate. Durable and chemical resistant body 5 cuff sizes are available for use (optional).



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265