

Product Name :
Sphygmomanometer Mercurial

Product Code :
PHARMACEUTICAL SUK-
AS03LAB058A



Description :

Sphygmomanometer Mercurial

Technical Specification :

Easy to read mercurial graduated column to 300mm Hg.

Mercurial blood pressure manometer, desk top type.

With coiled tubing approx. 152cm.

Velcro cuffs, adult and child.

Control valve and inflating bulb.

Stopcock for Hg tank.

Complete in aluminum case.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265