

Product Name :
Thermometer Laboratory

Product Code :
PHARMACEUTICAL SUK-
AS03LAB063A



Description :

Thermometer Laboratory

Technical Specification :

It is also a good idea to wrap the thermometer in a soft cloth to avoid breakage so you can safely conduct your science experiments. This thermometer fits in the holes of rubber stoppers. You might need to use a little cooking oil or Vaseline to lubricate the lab thermometers before they will fit.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265