

Product Name :
Double Unit Organ Bath With Thermostat

Product Code :
PHARMACEUTICAL SUK-
AS03LAB065A



Description :

Double Unit Organ Bath With Thermostat

Technical Specification :

This machine is provided with thermostatic temperature control, electric stirrer and related accessories for all kinds of research work. This product is manufactured using supreme grade raw material and advanced technology thus ensuring high quality and accuracy.

Technical Specifications :-

Replaceable metal sheathed heating element

Water Stirrer provided to evenly distribute the heat in the bath.

Sturdy metal frame

Bath: 12"x9"x8" with leak proof dove tailed joints

Digital Temperature Controller : Accuracy $\pm 1^{\circ}\text{C}$

Top of the bath band: 1"x3/8"

Supplied with all accessories made up of Brass (Non-Corrosive Metal).



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265