

**Product Name :**  
Stop Watch Digital Racer

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB071A



**Description :**

Stop Watch Digital Racer

**Technical Specification :**

The offered stop watch is designed in compliance with the set industry standards and assembled using components of superior quality to ensure high performance at user's end. Digital Stopwatch is a very versatile product. In normal mode it will display Time of the day in Hour, Minute, Seconds and the day of the week. Press Mode button once, and the stopwatch Mode will turn ON.

The product comes with a robust Plastic Storage Case. When readout becomes dim, it is time to replace the battery. Digital stop watch Racer is widely used for various applications in Sports, Laboratories and other industries. Stop watch racer is a handheld time piece designed to measure the amount of time elapsed from a particular time when it is activated to the time when the piece is deactivated.

**Features :-**

Large LCD Display for easy readout and Built-in Alarm Function.

Stopwatch : 1/100 second

Time Display : Hour, Minute, Seconds Calendar

Display : Day, Date, Month Alarm



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265