

Product Name :
Clamp Boss Heads

Product Code :
PHARMACEUTICAL SUK-
AS03LAB072A



Description :

Clamp Boss Heads

Technical Specification :

Metal thumbscrews with blue plastic heads allow for adjustment. For use with laboratory stands, this popular boss head is made of die casted alloy. It accepts rods up to 16 mm in diameter at a right angle to one another.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265