

**Product Name :**  
Projection Microscope

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB074A



**Description :**

Projection Microscope

**Technical Specification :**

Nosepiece Triple revolving nosepiece. Illumination 20 W Halogen bulb is provided with solid state variable light. Operates on 220 V AC 50 Hz. The structure of microscopic specimens visible to the naked eye by taking the projection on the screen during observation and demonstration to a group of students. Compact in design and provides high contrast and distortion free magnified images on fine ground screen.

**Specification :-**

Eyepieces Huygenian 10X and 15X.

Objectives Achromatic 5X, 10X, & 20 X (S). Screen 150 mm.

Diameters round glass projection screen rotatable through 360°.

Stage Fixed square stage size 115 x 115 mm with graduated detachable mechanical stage.

Accessories Cover, Duster, cleaning brush, spare bulb supplied in wooden case.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265