

Product Name :
Watch Glass

Product Code :
PHARMACEUTICAL SUK-
AS03LAB079A



Description :

Watch Glass

Technical Specification :

Watch glasses are made from flint glass and should not be heated directly. This 100 mm watch glass is useful as a beaker cover to prevent condensation escape and as an evaporating dish. Our large-size watch glasses have a slight depression, which makes them useful whenever a shallow vessel is required for a refraction surface.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265