

Product Name :
Electric Operated Centrifuges

Product Code :
PHARMACEUTICAL SUK-
AS03LAB080A



Description :

Electric Operated Centrifuges

Technical Specification :

It is being operated with the help of electricity where 220 volts current is applicable. Being reliable it is dual purpose equipment widely used in Gerber tests and solubility index tests.

Features :-

Precisely designed

Smooth functionality

Robustness



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265