

Product Name :
Sealing Machine

Product Code :
PHARMACEUTICAL SUK-
AS03LAB083A



Description :

Sealing Machine

Technical Specification :

The 8" Impulse Bag Sealer is the perfect one to Polyethylene, Mylar polyflex P.V.A polyvinylchloride and all other heat sealable materials. In addition, adjustable timer creates perfect seals for all kinds of bag thickness. You can adjust the time of the seal from the front knob, and when you press the top bar to seal, the red light will come on, and you will hear a click (this means the heating element is engaged). And as soon as the seal is done, the red light will go off. Ready to use instantly - No warm up needed.

Features :-

8" long sealing element.

Ready to use instantly - No warm up needed.

Can be used to seal kel-F, Mylar, polyflex, P.V.A, polyethylene, polyurethane, polyvinylchloride, polypropylene, pliofilm, tivac, saran.

Adjustable Sealing Control Timer with Ready Light With its adjustable timer, it will reach different plastic needs.

Labor-saving design with long handle.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265