

Product Name :
Bulk Density Apparatus

Product Code :
PHARMACEUTICAL SUK-
AS03LAB084A



Description :

Bulk Density Apparatus

Technical Specification :

The system has 'SET', 'UP' and 'DOWN' keys, for setting the total number of required strokes (Set Count) for the current measurement session. The system is equipped with microcontroller based digital display for displaying accurate and exact stroke count and soft touch membrane type keys for ease of operation.

User can set to any number of required strokes within the range from 1 to 9999 counts. The system automatically stores the last 'Set Count' value which is retained in the memory even when the instrument is switched off. A common lock cushioned holder for both the cylinders drops down smoothly from the height and holds the cylinders. Each stroke count is displayed on the 4 digit seven segment bright red LED display on the front panel.

Features :-

Highly Accurate and Easy to Operate

Microcontroller Based Digital Display

Extremely useful for Pharmaceutical Lab

Soft Touch Membrane Keys

3 LED's for status Indication

Audible Alarm Sound

Last Count Retention on Power Failure

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific