

**Product Name :**  
Antibiotic Zone Reader

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB086A



**Description :**

Antibiotic Zone Reader

**Technical Specification :**

The light from the inside source of the base passes up through the transparent & semitransparent of the agar then to a glass prism mounted at the magnified image of zone of inhibitions is clearly visible on the prism. The unit is measuring the dia meter of finished zone in a petri dish, for determining the strength of antibiotic material. This unit measures the dia meters of the inhibited zone within range of 0-35 mm with an accuracy of 0.2 mm.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

