

Product Name :
Mammalian Heart Perfusion Assembly

Product Code :
PHARMACEUTICAL SUK-
AS03LAB087A



Description :

Mammalian Heart Perfusion Assembly

Technical Specification :

The stirrer unit is completely independent, and dips into the bath in one corner, and fixes on its special stainless steel upright. A versatile assembly serving both for heart perfusion and as all-purpose thermostatic organ bath. Carrying the uprights and temperature of the bath is regulated and remains constant. Mammalian Heart Perfusion Assembly has a perspex connection box which is fitted to the bath frame, and carries all the electronic leads and connections.

Features :-

Strong build quality

Long lasting service life

Excellent thermostatic sensitivity



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265