

Product Name :
Bottle Sealing Machine

Product Code :
PHARMACEUTICAL SUK-
AS03LAB089A



Description :

Bottle Sealing Machine

Technical Specification :

Specifications :-

Conveyor Height: 750 mm to 840 mm Adjustable.

Container Size: 22 mm to 75 mm, Maximum Height 210 mm.

Output/Hour: 2400 to 4600 Containers (Depending on bottle size).

Electrical Motor: 1 HP/440 Volts/3 Phase/50 Hz (four wire) For Main Machine.

Cap Size: 22 mm to 38 mm (With the help of required change parts).

Overall Dimension: 1820 mm(L) X 750 (W) Approx.

Features :-

Built in front control panel

The unit is made compact & Versatile

S.S. Elegantly Matt Finished Body

S.S. Slat conveyor

Vibratory Bowl



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265