

Product Name :
Glassware Dryer

Product Code :
PHARMACEUTICAL SUK-
AS03LAB093A



Description :

Glassware Dryer

Technical Specification :

A time proportional high/low heat controlled with an automatic thermostatic safety cut off. It is used for drying laboratory glassware quickly and safety. By providing continuous warm air circulation through 12/24 hollow drying pegs with large or small bores with dust free filters labware is never too hot to handle.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265