

Product Name :
Tissue Homogenizer

Product Code :
PHARMACEUTICAL SUK-
AS03LAB094A



Description :

Tissue Homogenizer

Technical Specification :

It is lightweight enough to be hand-held and can also be mounted to a motor stand for lengthier processing. Tissue Homogenizer designed to handle a wide range of laboratory samples from 0.2 mL to 100 mL in volume. The usually process a tissue sample in less than 30 seconds and is compatible with a wide variety of generator probes for flexibility in application requirements.

The will accommodate a wide range of lab applications such as tissue homogenization, emulsions and suspensions, cell disruption, protein extraction and medical research.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265