

Product Name :
Piston Recorder

Product Code :
PHARMACEUTICAL SUK-
AS03LAB101A



Description :

Piston Recorder

Technical Specification :

It is examined by our quality associates on different characteristics to maintain its authenticity. The provided recorder is widely used in pharmaceutical industry to do practicals and researches in laboratories. Offered recorder is designed and manufactured with high precision in adherence with set market standards using optimum quality metal alloy and the latest technology.

Features :-

- Negligible maintenance
- Durable finish standards
- Accurate dimensions
- Sturdy construction



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265