

**Product Name :**  
Glass Bead Sterilizer

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB102A



**Description :**

Glass Bead Sterilizer

**Technical Specification :**

The unit is incorporated with a digital temperature indicator which continuously display the temperature being maintained in the glass beads. Glass Bead Sterilizer is known for its durable and sturdy construction. Manufactured using quality raw materials and modern machinery, this product is highly demanded among dental/ healthcare community. This product is used in hospitals, clinics, and several other such like establishments for the purpose of sterilizing various dental instruments, files, surgical instruments etc.

**Features :-**

- Bright Digital display for easy viewing and operations
- Ready to use within 60 seconds from cold start signalled with buzzer
- Latest Technology with highly accurate digital controller
- Low power consumption and Eco friendly



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265