

**Product Name :**  
Vibrating Reed

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB104A



**Description :**

Vibrating Reed

**Technical Specification :**

This unweighted ends are attached to a common base plate. The reed and base plate assembly is mounted in the field of an electromagnetic coil. This type of meter consists of a series of thin metal reeds of equal length usually of 4 mm wide and 0.5 mm thick. But slightly different weights are attached to one end of the reeds. This added weights made each reed mechanically resonant at a particular frequency. The weighted ends of the reeds are visible through a narrow window in the face of the meter housing and are painted white.

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com