

**Product Name :**  
Suppository Mould

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB106A



**Description :**

Suppository Mould

**Technical Specification :**

Nickel plated various capacities with a number of moulds in each set.  
Suppository Mould Made from brass/gun metal machine made perfect shapes  
smooth and polished cavities.

**Capacity :-**

- 1 gm, 1.5 gms & 2 gms (15.4 grains)
- 4-moulds set
- 6-moulds set
- 12 moulds set.



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265