

Product Name :
Student Spirometer

Product Code :
PHARMACEUTICAL SUK-
AS03LAB109A



Description :

Student Spirometer

Technical Specification :

The chamber is made with stainless steel corrugated rubber tube and mouth piece. The product has a chain attached which compensates for counter balance to flat through pulley calibrated to denote volume inlet and outlet tubes.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265