

**Product Name :**  
Respiration Pump

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB116A



**Description :**

Respiration Pump

**Technical Specification :**

The stroke is adjustable and is graduated up to 500 cc. Motor fixed on adjustable castor wheeled base. Single cylinder pump with mechanically operated valve. All metal construction. The driving cone has 3 slots of 6", 9", 12" diameter and is accurately balanced. A 0.25 hp motor workable on 220v ac 50hz single phase operates with the pump.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265