

Product Name :
Perimeter Priestly Smith

Product Code :
PHARMACEUTICAL SUK-
AS03LAB117A



Description :

Perimeter Priestly Smith

Technical Specification :

A metallic graduated scale with a movable pin-punch pointer is fixed to mark on the perimeter chart paper. The perimeter is an instrument for accurate charting of the field of vision. It consists of. A vertical stand on which a metallic arc is pivoted. It also bear a circular black disc to read the meridian in which the arc in shape of a semicircle with radius 330 mm. Its rotating in any direction and fix at any position with a tightening screw.

The arc is graduated from 0° to 90° with a movable test object which is a white color spot of size 10mm diameter painted on a circular disc it can be moved with arc. An adjustable chin rest and a detachable leveling bar is fixed in front of the metallic arc. A scale and circular chart frame to hold the chart paper in position is provided at the back of the circular disc.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265