

Product Name :
Magnetic Stirrer With Hot Plate

Product Code :
PHARMACEUTICAL SUK-
AS03LAB119A



Description :

Magnetic Stirrer With Hot Plate

Technical Specification :

Stirring and heating units are properly insulated. In model with heating system temp. is controlled by an energy regulator. The provided stirrer is a laboratory device that employs a rotating magnetic field to cause a stir bar immersed in a liquid to spin very quickly, thus stirring it. Compact stirring device utilizing a rotating magnetic field applied from below the container. The stirring is accomplished by a teflon rotor placed in the liquid to be stirred. Complete with pilot lamp, plug cord, continuously variable speed regulator. The body is made of mild steel duly finished in stoving enamel and mounted with a stainless steel top which serves as a hot plate. An electronic speed controller is also provided. The equipment works on 220V AC 50 HZ single phase. This stirrer is used in schools, colleges, research centers etc.

Features :-

Superior efficiency

Low maintenance

User friendly design

Maintenance free



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265