

**Product Name :**  
Laminar Air Flow

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB120A



**Description :**

Laminar Air Flow

**Technical Specification :**

The Laminar Air Flow also has an acrylic panel, fluorescent tube light & pressure differential manometer inside shelf that is made of stainless steel. A Laminar Air Flow is designed to provide working space comprising as standard features. The Laminar Air Flow is a HEPA filter having efficiency of 99.97% down to 0.3% microns, prefilter 90% down to 5 microns, suitable motor blowing assemblies, working table top side.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265